

# Mayer's Hematoxylin Solution

*ready-to-use*

**42**  
life sciences

*WE HAVE THE SOLUTION!*

*Ready-to-use Mayer's Hematoxylin Solution in stabilizing buffer:*

- one stable component*
- optimized for immunohistochemical (IHC) procedures*
- non-toxic*




## Introduction

Mayer's Hematoxylin Solution is a staining reagent for histological applications. Nuclei in stained sections will be a blue-violet color with crisp nuclear detail. This reagent is suitable for use with all chromogens commonly used in immunohistochemical (IHC) applications.

## Description

Ready-to-use Mayer's Hematoxylin Solution in stabilizing buffer. Does not contain chloralhydrate. Optimized for the use in immunohistochemical (IHC) procedures.



Counterstain with *Mayer's Hematoxylin Solution* on tonsil section.

## *Order Information*

Prod. No.	Product
S-99025-20	Mayer's Hematoxylin Solution - 20 ml
S-99025-1000	Mayer's Hematoxylin Solution - 1000 ml

42 life sciences GmbH & Co. KG

Fischkai 1  
D-27572 Bremerhaven  
Germany  
[www.42ls.com](http://www.42ls.com)  
[info@42ls.com](mailto:info@42ls.com)